

FLUORESCENCE-BIOSENS Line Fluorescence BioSensor Systems



System of Measurement for Algae (SMA) is a small portable instrument to optically characterize the fluorescence emission and the optical density (correlated to a cell density) of algae culture (e.g. *Haematococcus pluvialis*, *Spirulina sp.*, *Chlorella minutissima*, etc.). It is possible to connect the SMA instrument directly to a by-pass of a bioreactor to monitor the algae culture condition.

Keyboard and display can be used to set-up the parameters.

It is important to verify the cleaning procedures above all for OD measurements and the dark time for the fluorescence emission.




The biological material is introduced into a Delrin tube by a solenoid injection valve; optical sensors allow to measure the fluorescence emitted by chlorophyll inside the green algae cells in culture medium and the density of these cells culture medium through a light transmission measurement. It is equipped with current measurement system, and an electronic control board for data read-out, processing and PC transfer.

Main Characteristics

Functional Characteristic	Optical density and fluorescence biosensor used in environmental analysis
Sample	Algae culture monitoring
Readings	This instrument measures the Optical Density (O.D.) detected by the optical system and the fluorescence emission
Optical	
Fluorescence Cell	Light source LED, peak emission wavelength: 475 nm
Measurement Optical Density (Cells density measurement)	Light source LED emission peak wavelength: 750 nm
Instrument Dimensions	22 x 15 x 10 cm (W x L x H)
Weight	980 gr

Fluorescence Test on green algae (*Haematococcus pluvialis*) Dark time 10 min; Light source 10 sec; Level led 127; Relax time 2min



It is possible to customize the instrument meet the customer requirement (e.g. Leds wavelength, Leds intensity, etc).

CONTACT:

Biosensor s.r.l.
Via degli Olmetti 44, 00060 Formello (Rome) –
Italy

tel/fax: +39 06 9075116
e-mail: info@biosensor.it
website: www.biosensor.it